



## Datasheet

### Panther 5 Series

Revision B



## Technical Specifications

### General Information

|                           |   |
|---------------------------|---|
| Intended Use              | The PANTHER 5 ventilator is: <ul style="list-style-type: none"> <li>Intended for respiratory treatment in invasive and non-invasive of neonatal<sup>(N)</sup>, pediatric and adult patients</li> <li>Used in hospitals, professional healthcare facilities and transportation of patients within such facilities</li> </ul> |
| Instructions for use      | Please read the Panther 5 operator's manual   |
| Legal Manufacturer        | Origin Medical Devices Inc.   |
| Size (W x D x H)          | 350 x 360 x 390 mm<br>13.7" x 14.2" x 15.4"   |
| Weight (Ventilator)       | 11 Kg (24 lbs)  |
| Power                     | 110 to 240 V AC<br>50 to 60 Hz  |
| Internal Battery          | Li-ion, 98Wh, 14.4V   |
| Operating Time on Battery | > 3 hours under standard conditions   |
| Recharge Time             | Approximately 3 hours   |

### Oxygen Supply and Monitoring

|                     |  |
|---------------------|--|
| High Pressure Range | 35 to 87 PSI   |
| Connector Type      | DISS 1240, NIST or other per region  |
| Low Pressure        | Low flow/pressure inlet  |
| Monitoring          | O <sub>2</sub> sensor on outlet.<br>Galvanic or Paramagnetic options available |

<sup>(N)</sup> feature is currently not available in the USA

### Operational

|                       |  |
|-----------------------|--|
| Enclosure Rating      | IP22                                   |
| Operating Temperature | 10 to +40°C                            |
| Operating Humidity    | 10 to 90% Non-Condensing               |
| Storage Temperature   | -20 to +60°C                           |
| Storage Humidity      | 10 to 90% Non-Condensing               |
| Barometric Pressure   | 700 to 1060 kPa internally compensated |
| Altitude - operation  | 0 to 3,280 m<br>(0 to 10,000 ft)       |

### Flow and Pressure

|                |                                |
|----------------|--------------------------------|
| Pressure Range | -50 to +100 cmH <sub>2</sub> O |
| Flow Range     | 0 to 240 lpm                   |

### Functionality and safety standards



Complies with requirements and classification IIb of Medical Device Directive 93/42/EEC.

ISO 80601-2-12:2011  
ISO 60601-1-2:2014  
EN 60601-1  
ISO 60601-1-8:2007 + A11:2017  
IEC 60601-2-49:2011  
ISO 80601-2-55:2018

## User Interface

|                              |   |
|------------------------------|---|
| Display                      | 15" TFT with PCAP touchscreen                                       |
| Control Interface            | Touchscreen<br>Encoder knob with LED                                |
| Audible Indicators           | Speaker and Buzzer  |
| Additional Visual Indicators | RED, YELLOW, GREEN indicators for alarms, ventilation in power save |
| Additional Visual Sensors    | Ambient light detector for automated display intensity control      |

## Mode Selections

|       |  |
|-------|--|
| Types | Invasive<br>Non-invasive<br>High flow O <sub>2</sub> therapy |
|-------|--|

## Invasive Ventilation Modes

|  |  |
|--|--|
| Controlled Ventilation   | Volume Control (VC)<br>Pressure Control (PC)<br>Pressure regulated volume flowcontrol (PRVC) |
| Support Ventilation  | Pressure Support (PS)<br>Volume Support (VS)   |
| SIMV (synchronized intermittent mandatory ventilation)   | VC + VS<br>VC + PS<br>PC + PS<br>PRVC + PS<br>PRVC + VS                                      |
| SMART Mode™ (Automatic transition between a control mode and a spont mode based on presence or absence of patient efforts) | VC + VS<br>VC + PS<br>PC + PS<br>PRVC + PS<br>PRVC + VS                                      |
| BiLevel  | Dual PEEP with two defined PS levels   |
| APRV   | Dual PEEP CPAP   |

## Non-Invasive Ventilation Modes

|  |  |
|--|--|
| Controlled Ventilation                                 | Volume Control (VC)<br>Pressure Control (PC)                                       |
| Support Ventilation                                    | Pressure Support (PS)  |
| SIMV (synchronized intermittent mandatory ventilation) | VC + PS<br>PC + PS   |
| SMART Mode™  | VC + PS<br>PC + PS   |
| NIPPV <sup>(N)</sup>                                   | Nasal Intermittent Positive Pressure Ventilation (synchronized to patient efforts) |

## Non-Invasive compensation

|   |   |
|---|---|
| Non-invasive max leak compensation        | Adult <ul style="list-style-type: none"> <li>Inspiratory: 200 lpm</li> <li>Expiratory: 60 lpm</li> </ul> Pediatric and neonatal <ul style="list-style-type: none"> <li>Expiratory 25 lpm</li> </ul> |
| Inspiratory volume compensation during VC | User selectable: ON/OFF<br>When ON, volume is compensated for up to twice the defined volume  |

## Invasive Compensation

|   |  |
|---|--|
| Invasive max leak compensation when activated | Adult <ul style="list-style-type: none"> <li>Inspiratory: 200 lpm</li> <li>Expiratory: 25 lpm</li> </ul> Pediatric and neonatal<br>Expiratory 25 lpm |
| Inspiratory volume compensation during VC     | User selectable: ON/OFF<br>When ON, volume is compensated for up to an additional 50% of defined volume  |

## High Flow O<sub>2</sub> therapy

|        |  |
|--------|--|
| Oxygen | 21 to 100%   |
| Flow   | Adult: 1 to 80 lpm<br>Pediatric: 1 to 80 lpm<br>Neonate: 1 to 25 lpm |

## SBT (spontaneous breathing trial)

|                           |   |
|---------------------------|---|
| SBT Time                  | 15 to 120 Minutes   |
| Oxygen (O <sub>2</sub> %) | 21 to 100%  |
| PEEP                      | 0 to 40 cmH <sub>2</sub> O  |
| Support Pressure          | 0 to 75 cmH <sub>2</sub> O  |
| Support Slope             | 1 to 10   |
| Termination               | Analyzes patient degradation using a variety of indications and automatically resume normal ventilation when needed |

## Additional Functions

|  |   |
|--|---|
| Speaking Valve                                   | Automatically adjusts specific alarms, disables activation of conflicting features to enable safe use of a speaking valve. When turned off reactivates disabled features and returns alarms to normal |
| Demand Flow (in VC only)                         | Detects patient's need of additional flow and automatically transitions to PS for that specific breath  |
| Auto exhalation Sensitivity (E <sub>SENS</sub> ) | Automatic Breath-by-Breath ventilator management of the exhalation sensitivity setting  |
| SMART Trigger™                                   | Proprietary triggering mechanism which significantly improves trigger detection in high and varying leaks as well as COPD patients  |

## IBW Calculation

|               |   |
|---------------|---|
| Weight Ranges | Adult: 25 to 144 Kg<br>Pediatric: 2.9 to 24 Kg<br>Neonate: 0.4 to 2.8 Kg  |
| High Ranges   | Adult: 125 to 256 cm<br>Pediatric: 48 to 124 cm<br>Neonate: 26.5 to 47 cm |
| Units         | Allows entry of height in cm or inches and weight in Kg or lbs            |
| Gender        | Male or female  |

## Display Configurations

|                 |   |
|-----------------|---|
| Waveforms       | <ul style="list-style-type: none"> <li>• Circuit Pressure</li> <li>• Flow</li> <li>• Volume</li> <li>• Circuit + Tracheal Pressure (TC)</li> <li>• CO<sub>2</sub></li> <li>• SpO<sub>2</sub></li> </ul>   |
| Loops           | <ul style="list-style-type: none"> <li>• Pressure Volume (PV)</li> <li>• Flow Volume (FV)</li> <li>• Volume Flow (VF)</li> <li>• Single Breath CO<sub>2</sub> Curve</li> </ul>  |
| Reference Loops | Shows up to two out of eight saved loops superimposed on the live loops along with event information prior to the loop save   |
| Trends          | <ul style="list-style-type: none"> <li>• Shows one live waveform, two selected trends and 25 monitored values corresponding to the trend cursor position.</li> <li>• Trend views can be selected from 25 trended parameters which are recorded per breath (no faster than once a second)</li> <li>• Views can be zoomed and scrolled with the x-axis or finger swipe</li> <li>• Trends record 72 hours of data</li> </ul> |

## Maneuvers

### P0.1 (P100)

User initiated automated maneuver to measure the patient's respiratory drive during the first 100 ms of inspiratory effort when the airway is occluded.

|          |  |
|----------|--|
| Max Time | Adult: 8 seconds<br>Pediatric: 6 seconds |
|----------|--|

### Smart NIF

User initiated automated maneuver to measure the patient's inspiratory muscle strength during airway occlusion. Provides visual and audible indications and automated analysis of patient fatigue to increase safety and patient synchronization.

|                   |  |
|-------------------|--|
| Max maneuver Time | Adult: 20 seconds<br>Ped/Neo: 10 seconds |
|-------------------|--|

|            |  |
|------------|--|
| Initiation | Audible BEEP and LED flashing indicated to the patient and clinician |
|------------|--|

|             |   |
|-------------|---|
| Termination | Time and automated detection of patient fatigue |
|-------------|---|

### Suction Maneuver

Provides an automated safe management of suction procedure by automatically transitions between suction phases, adjusts alarms, and detects reconnection to resume ventilation.

|       |  |
|-------|--|
| Types | <ul style="list-style-type: none"> <li>• Open Suction</li> <li>• Closed Suction</li> </ul> |
|-------|--|

|                         |           |
|-------------------------|-----------|
| O <sub>2</sub> settings | 21 – 100% |
|-------------------------|-----------|

|        |   |
|--------|---|
| Phases | <ul style="list-style-type: none"> <li>• Pre-Oxygenation at set O<sub>2</sub> level</li> <li>• Suction: auto detected on open suction, manual press in closed suction</li> <li>• Post-Oxygenation: auto detected on open suction, manual press in closed suction</li> <li>• Automated phase timeouts</li> </ul> |
|--------|---|

|         |  |
|---------|--|
| Effects | Oxygen, alarms and delayed activation of alarms automated by the system. |
|---------|--|

### PV (Slow Inflation/Deflation)

A slow inflation/deflation PV Maneuver is both a diagnostic and therapeutic tool that provides information that may be used to optimize PEEP, tidal volume and other ventilator settings to allow lung protective ventilation. Upon maneuver completion, ventilation transitions back to the settings prior to the maneuver at the user defined End PEEP setting.

|            |                            |
|------------|----------------------------|
| Start PEEP | 0 to 40 cmH <sub>2</sub> O |
|------------|----------------------------|

|              |                 |
|--------------|-----------------|
| PEEP EQ Time | 0.0 to 30.0 sec |
|--------------|-----------------|

|                     |                                 |
|---------------------|---------------------------------|
| Inflation/Deflation | 2 to 5 cmH <sub>2</sub> O / sec |
|---------------------|---------------------------------|

|                 |                            |
|-----------------|----------------------------|
| Target Pressure | 5 to 60 cmH <sub>2</sub> O |
|-----------------|----------------------------|

|                 |                 |
|-----------------|-----------------|
| Pause at Target | 0.0 to 30.0 sec |
|-----------------|-----------------|

|          |                            |
|----------|----------------------------|
| End PEEP | 0 to 40 cmH <sub>2</sub> O |
|----------|----------------------------|

|            |        |
|------------|--------|
| Time limit | 60 sec |
|------------|--------|

|                    |  |
|--------------------|--|
| Safety Termination | Time and on patient effort (resumes ventilation) |
|--------------------|--|

|       |  |
|-------|--|
| Views | <ul style="list-style-type: none"> <li>• PV loop during the maneuver, auto scaled</li> <li>• Maneuver graph shows expected maneuver progress and progress during the maneuver</li> </ul> |
|-------|--|

|             |   |
|-------------|---|
| Measurement | Four cursors for four inflection points |
|-------------|---|

### Recruitment

Single or multi step recruitment maneuver (RM) via continuous ventilation at user defined step settings. Upon maneuver completion, ventilation transitions back to the settings prior to the maneuver at the user defined End PEEP setting.

|                 |         |
|-----------------|---------|
| Number of steps | 1 to 20 |
|-----------------|---------|

|        |                 |
|--------|-----------------|
| T High | 1.0 to 59.0 sec |
|--------|-----------------|

|       |                |
|-------|----------------|
| T Low | 1.0 to 5.0 sec |
|-------|----------------|

|        |                             |
|--------|-----------------------------|
| P High | 10 to 40 cmH <sub>2</sub> O |
|--------|-----------------------------|

|       |                            |
|-------|----------------------------|
| P Low | 0 to 30 cmH <sub>2</sub> O |
|-------|----------------------------|

|          |                            |
|----------|----------------------------|
| End PEEP | 0 to 30 cmH <sub>2</sub> O |
|----------|----------------------------|

|       |   |
|-------|---|
| Views | Graphical representation of the maneuver and its progress |
|-------|---|

## Direct Access Functions

|                         |  |
|-------------------------|--|
| Elevated O <sub>2</sub> | User adjustable O <sub>2</sub> level active for up to 120 seconds      |
| Manual Breath           | Activates a mandatory breath upon pressing during the expiratory phase |

## Tube Compensation

|            |  |
|------------|--|
| Activation | ON/OFF (restrictions apply)  |
| Tube Type  | <ul style="list-style-type: none"> <li>• Endotracheal</li> <li>• Tracheostomy</li> </ul> |
| Tube ID    | Adult: 5.5 to 10.0 mm<br>Pediatric: 4.0 to 6.5 mm<br>Neonate: 2.0 to 4.5 mm              |
| Length     | Adult: 2.0 to 30.0 cm<br>Pediatric: 2.0 to 26.0 cm<br>Neonate: 2.0 to 15.0 cm            |
| Support %  | 10 to 100%   |

## Nebulization

### Pneumatic

|                |  |
|----------------|--|
| Flow           | 7 lpm, Oxygen                                    |
| Operating Time | 5, 10, 20, 30 minutes                            |
| Compensation   | Volume is compensated for the added flow         |
| Automation     | Automated termination under violating conditions |

### Aerogen

|                 |  |
|-----------------|--|
| Method          | Direct drive to nebulizer  |
| Supported Types | SOLO or PRO  |
| Controls        | <ul style="list-style-type: none"> <li>• Selection of type</li> <li>• Continuous option</li> <li>• Run time and Extend time</li> </ul> |
| Visuals         | <ul style="list-style-type: none"> <li>• Time since start</li> <li>• Time to run</li> <li>• Operational status</li> </ul>              |

## Capnography

|              |   |
|--------------|---|
| Measurements | Measures EtCO <sub>2</sub> and real time inspired and expired CO <sub>2</sub> |
|--------------|---|

### Calculated Parameters

|  |   |
|--|---|
| End-Tidal CO <sub>2</sub>                                | EtCO <sub>2</sub>   |
| Fractional End-Tidal CO <sub>2</sub> Concentration       | FetCO <sub>2</sub>  |
| Partial Pressure of Mean Expired CO <sub>2</sub>         | PeCO <sub>2</sub>   |
| Fractional Concentration of Mean Expired CO <sub>2</sub> | FeCO <sub>2</sub>   |
| Exhaled CO <sub>2</sub> Volume                           | V <sub>TE</sub> CO <sub>2</sub>   |
| Inspired CO <sub>2</sub> Volume                          | V <sub>TI</sub> CO <sub>2</sub>   |
| Exhaled Volume of CO <sub>2</sub> Per Minute             | VCO <sub>2</sub>  |
| Alveolar Tidal Volume                                    | Valv  |
| Alveolar Minute Volume                                   | Valv/min  |
| Anatomical Dead Space                                    | V <sub>d</sub> ana  |
| Alveolar Dead Space                                      | V <sub>d</sub> alv  |
| Physiological Dead Space to Tidal Volume Ratio           | VD/VT Eng   |
| Physiological Dead Space to Tidal Volume Ratio           | VD/VT Bohr Est  |
| Graphical  | <ul style="list-style-type: none"> <li>• Single breath CO<sub>2</sub> curve</li> <li>• Realtime exhaled CO<sub>2</sub> over time</li> </ul> |

## Oximetry

|              |   |
|--------------|---|
| Measurements | <ul style="list-style-type: none"> <li>• SpO<sub>2</sub></li> <li>• Heart Rate</li> <li>• SpO<sub>2</sub>/O<sub>2</sub> (Ratio Approximation to PaO<sub>2</sub>/FiO<sub>2</sub>)</li> <li>• Signal Level</li> </ul> |
|--------------|---|

## Log

|                    |   |
|--------------------|---|
| Logged Information | <ul style="list-style-type: none"> <li>• Changes</li> <li>• Alerts</li> <li>• Operations</li> </ul> |
| Number of Entries  | 5,000   |

## Settings

|   |  |
|---|--|
| PEEP / CPAP                               | 0 to 40 cmH <sub>2</sub> O   |
| Pressure (PC)<br>(above PEEP)             | 5 to 90 cmH <sub>2</sub> O   |
| Support Pressure<br>(PS) above PEEP       | 0 to 75 cmH <sub>2</sub> O   |
| Mandatory Slope<br>(Mand Slope)           | 1 to 10<br>(1 is the fastest)  |
| Spontaneous Slope<br>(Support Slope)      | 1 to 10<br>(1 is the fastest)  |
| Exhalation<br>Sensitivity (Esens)         | 5 to 80 %  |
| Max Spont Breath<br>Time (Support TI)     | Adult: 0.4 to 5.0 sec<br>Pediatric: 0.4 to 3.0 sec<br>Neonate: 0.2 to 2.0 sec                                  |
| Tidal Volume (VT)<br>Range                | Adult: 100 to 2500 ml<br>Pediatric: 20 to 500 ml<br>Neonate <sup>(N)</sup> : 2 <sup>(1)</sup> to 100 ml        |
| Waveform (VC)                             | Square / Decelerating  |
| Inspiratory Time (TI)                     | Adult: 0.10 to 5.00 sec<br>Pediatric: 0.10 to 4.00 sec<br>Neonate: 0.10 to 3.00 sec<br>BiLevel: Up to 59.8 sec |
| Plateau Time<br>(Insp Hold)               | 0.0 to 3.0 sec   |
| Respiratory Rate                          | Adult: 1 to 110 b/min<br>Pediatric: 1 to 120 b/min<br>Neonate: 1 to 150 b/min                                  |
| SMART Time<br>(Applies to SMART<br>Mode™) | 3 to 15 sec  |
| Rate NIPPV                                | 1 to 150 b/min   |
| Pressure Trigger                          | -15 to -0.1 cmH <sub>2</sub> O   |
| Flow Trigger                              | 0.1 to 20 lpm  |
| SMART Trigger™                            | 1 to 7   |
| Oxygen (O <sub>2</sub> %)                 | 21 to 100%   |

<sup>(1)</sup> 5ml in VC, 2ml in pressure modes

<sup>(N)</sup> feature is currently not available in the USA

## Apnea

|                               |  |
|-------------------------------|--|
| Apnea Time                    | 0 to 60 sec and OFF<br>(OFF allowed in SPONT<br>when PS is set ≤ 5 cmH <sub>2</sub> O) |
| Oxygen (O <sub>2</sub> %)     | 21 to 100%   |
| Inspiratory Time (TI)         | Adult: 0.10 to 5.00 sec<br>Pediatric: 0.10 to 4.00 sec<br>Neonate: 0.10 to 3.00 sec    |
| Pressure (PC)<br>(above PEEP) | 5 to 90 cmH <sub>2</sub> O   |
| Respiratory Rate              | Adult: 1 to 110 b/min<br>Pediatric: 1 to 120 b/min<br>Neonate: 1 to 150 b/min          |
| Slope                         | Uses set Slope   |
| Triggers                      | Uses set triggers  |

## Operating Ranges

|           |         |
|-----------|---------|
| Peak Flow | 240 lpm |
|-----------|---------|

## Monitored Parameters

### Pressure

|                                    |                         |
|------------------------------------|-------------------------|
| Peak Pressure                      | P <sub>peak</sub>       |
| End Inspiratory Pressure           | P <sub>Insp</sub>       |
| Mean Airway Pressure               | P <sub>mean</sub>       |
| End Expiratory Pressure            | PEEP                    |
| Calculated Tracheal Pressure       | P <sub>trach Insp</sub> |
| Plateau Pressure                   | P <sub>Plateau</sub>    |
| Intrinsic (Auto) PEEP              | Intr. PEEP              |
| Total PEEP                         | PEEP Tot                |
| Delta Air Pressure                 | dPAW                    |
| Breathing Drive Occlusion Pressure | P <sub>O.1</sub>        |
| Negative Inspiratory Pressure      | NIF                     |

### Volume

|                                       |                       |
|---------------------------------------|-----------------------|
| Inspired Tidal Volume                 | V <sub>Ti</sub>       |
| Exhaled Tidal Volume                  | V <sub>TE</sub>       |
| Spontaneous V <sub>Ti</sub>           | Spont V <sub>Ti</sub> |
| Spontaneous V <sub>TE</sub>           | Spont V <sub>TE</sub> |
| VTI Normalized to Patient Body Weight | V <sub>Ti</sub> /PBW  |
| VTE Normalized to Patient Body Weight | V <sub>TE</sub> /PBW  |
| Inspired Minute Volume                | Ṁ <sub>i</sub>        |
| Exhaled Minute Volume                 | Ṁ <sub>e</sub>        |
| Spont Inspired Minute Volume          | Spont Ṁ <sub>i</sub>  |
| Spont Exhaled Minute Volume           | Spont Ṁ <sub>e</sub>  |

### Flow and Leak

|                         |                      |
|-------------------------|----------------------|
| Peak Inspiratory Flow   | PIF                  |
| Peak Expiratory Flow    | PEF                  |
| Delivered Oxygen        | O <sub>2</sub>       |
| Inspiratory Leak (lpm)  | Insp Leak lpm        |
| Inspiratory Leak (%)    | Insp Leak %          |
| Average Total Leak Rate | Avg Leak lpm         |
| Inspiratory Leak Volume | V <sub>leak</sub> ml |

### Rate and Timing

|                                    |          |
|------------------------------------|----------|
| Total Breath Rate                  | Total BR |
| Mandatory Respiratory Rate         | Mand BR  |
| Spontaneous Respiratory Rate       | Spont BR |
| Inspiratory Time                   | Last Ti  |
| Expiratory Time                    | Te       |
| Spontaneous Inspiratory Time Ratio | TH/Ttot  |
| Ratio between THigh and TLow       | TH:TL    |
| Inspiratory to Expiratory Ratio    | I:E      |

### Mechanics

|                                   |                    |
|-----------------------------------|--------------------|
| Inspiratory Pressure Time Product | PTP                |
| Static Compliance                 | C <sub>stat</sub>  |
| Dynamic Compliance                | C <sub>dyn</sub>   |
| Static Resistance                 | R <sub>Stat</sub>  |
| Expiratory Resistance             | RE                 |
| Inspiratory Time Constant         | RC <sub>Insp</sub> |
| Expiratory Time Constant          | RC <sub>Exp</sub>  |
| Rapid Shallow Breathing Index     | RSBI               |
| Work of Breathing Imposed         | WOB Imposed        |

### Capnography and Oximetry

See Capnography and Oximetry sections above.



## Adjustable Alarms

|                                       |   |
|---------------------------------------|---|
| Pressure High                         | 6 to 100 cmH <sub>2</sub> O   |
| Pressure Low                          | 3 to 97 cmH <sub>2</sub> O  |
| Minute Volume (V̇ <sub>e</sub> ) High | 0.5 to 100 lpm Adult<br>0.5 to 30 lpm Pediatric<br>0.5 to 10 lpm Neonate                            |
| Minute Volume (V̇ <sub>e</sub> ) Low  | OFF to 0.1 to 99.5 lpm Adult<br>OFF to 0.05 to 29.5 lpm Pediatric<br>OFF to 0.01 to 9.5 lpm Neonate |
| V <sub>TE</sub> High                  | 25 to 3000 to OFF ml Adult<br>25 to 700 to OFF ml Pediatric<br>5 to 300 to OFF ml Neonate           |
| V <sub>TE</sub> Low                   | OFF to 1 to 2500 ml Adult<br>OFF to 1 to 690 ml Pediatric<br>OFF to 1 to 295 ml Neonate             |
| Spont V <sub>TE</sub> High            | 25 to 3000 to OFF ml Adult<br>25 to 700 to OFF ml Pediatric<br>5 to 300 to OFF ml Neonate           |
| Spont V <sub>TE</sub> Low             | OFF to 1 to 2500 ml Adult<br>OFF to 1 to 690 ml Pediatric<br>OFF to 1 to 295 ml Neonate             |
| Rate High                             | 10 to 110 b/min Adult<br>10 to 130 b/min Pediatric<br>10 to 170 b/min Neonate                       |
| Rate Low                              | 1 to 109 b/min Adult<br>1 to 129 b/min Pediatric<br>1 to 169 b/min Neonate                          |
| Disconnect Sensitivity (Dsens)        | 20 to 95 %  |
| V <sub>Ti</sub> Limit                 | 105 to 3000 ml Adult<br>25 to 750 ml Pediatric<br>6 to 300 ml Neonate                               |

|                         |                       |
|-------------------------|-----------------------|
| ETCO <sub>2</sub> High  | 10 to 150 to OFF mmHg |
| ETCO <sub>2</sub> Low   | OFF to 5 to 60 mmHg   |
| VteCO <sub>2</sub> High | 0.2 to 100 to OFF ml  |
| VteCO <sub>2</sub> Low  | OFF to 0.1 to 99 ml   |
| SBT Rate High           | 5 to 80 to OFF b/min  |
| SBT Rate Low            | OFF to 1 to 75 b/min  |
| SBT RSBI High           | 5 to 900 to OFF       |
| SBT RSBI Low            | OFF to 5 to 895       |
| Leak High               | 5 to 95 %             |
| SPO <sub>2</sub> High   | 71 to 100%            |
| SPO <sub>2</sub> Low    | 70 to 99%             |
| Heart Rate High         | 45 to 245 bpm         |
| Heart Rate Low          | 40 to 240 bpm         |

## Non-Adjustable Alarms

|                                  |                      |
|----------------------------------|----------------------|
| Standby                          | Occlusion            |
| Low PEEP                         | High PEEP            |
| PRVC Limited by High P           | VS Limited by High P |
| Circuit Open                     | Apnea                |
| Low O <sub>2</sub>               | High O <sub>2</sub>  |
| No O <sub>2</sub> Inlet Pressure | Aerogen Fault        |
| Battery Gauge Error              | Battery Hot          |
| Battery Low                      | Battery Empty        |
| Shutting Down                    | Charger Fault        |

### Additional Technical Alarms

Additional CO<sub>2</sub> module related error alarms

Additions SpO<sub>2</sub> module related error alarms

## Communication Interfaces

### Serial RS232

- Sends automatic data to nurse call station
- Can be configured to send the required data under different conditions
- Software plug-ins for required protocols

### Ethernet

- Sending automatic data as well as provides online monitoring, log reading and remote control
- Software enables connection to dedicated control/monitoring software that run on remote computers/tablets/phones or standard control centers
- USB Host connection for saving of logs, screen images and uploading software updates from standard USB memory sticks

## External Interfaces

Capnography Module

SpO<sub>2</sub> Modules

Direct Aerogen Nebulizer

Dry contact remote alarm connections with/without cable disconnection detection

## Ventilator Options

### Software

|                                  |   |
|----------------------------------|---|
| Neonatal Suite <sup>(N)</sup>    | Neo patient type, NIPPV mode  |
| Synchrony Suite                  | Demand Flow<br>Auto E <sub>SENS</sub><br>Tube Comp                  |
| Diagnostic and Therapeutic Suite | PV Maneuver<br>Recruitment Maneuver<br>Suction Maneuver             |
| Extended Monitoring Suite        | Volumetric Capnography and Oximetry monitoring, alarms and trending |

### Hardware

|                       |                     |
|-----------------------|---------------------|
| O <sub>2</sub> Sensor | Paramagnetic sensor |
|-----------------------|---------------------|

The ventilator includes all hardware to fully support all features and all software options. There is no need to install additional internal hardware for any option.

<sup>(N)</sup> feature is currently not available in the USA